


# Optimising Vitrification

**Dave Morroll, Director of Clinical Support**

**david.morroll@coopersurgical.com**

---

Israel, February 2019



# Aims & Objectives

- introduction to cryobiology and use of CPAs
- highlight key areas of concern
- detail stage-by-stage considerations
- discussion

# Cryopreservation

- the main aim of any cryopreservation programme is to ensure a high percentage of survival and viability after thawing
- six steps:
  - equilibration with cryoprotectant agents to induce dehydration
  - reduction of the temperature to a cooling point
  - storage
  - warming
  - rehydration
  - resumption of cellular functions







# Cryoprotectants



- agents which influence the transition from water to ice
- stabilise bio-molecules as water substitutes
- reduce ice nucleation and growth during thawing
- intracellular or permeating
- extracellular or non-permeating

# Cryoprotectants

- intracellular or permeating





- enter cells and replace water
- as they enter cells more slowly than water exits, there is shrinking and re-expansion



# Cryoprotectants

- extracellular or non-permeating



# Warming/thawing

- after thawing, need progressive dilution of CPA
- form osmotic gradient by creating an extracellular environment with lower osmolarity
- this results in the gradual movement of water in to the cell = rehydration
- associated diffusion of intracellular CPA to the extracellular environment

# Optimising vitrification: what is important?



# Recurring questions...

- what should be vitrified?
- which media/CPA is best?
- how and when must we control temperature?
- are volumes of media important?
- how long must we equilibrate?
- must we collapse blastocysts and, if so, how is equilibration affected?
- choice of carrier: open or closed systems?
- how do we optimise warming procedures?

# Embryo grading: what to vitrify?

- cleavage stage
  - appropriate cell number (by age hpi)
  - <20% blast diameter difference; <20% fragmentation
- blastocyst stage
  - full blastocyst
    - compacted ICM, loosely adhered, some individual cells
    - confluent TE layer, fewer and larger cells
    - **leave early blasts to D6 but FBT as a day 5**

# Embryo grading


## Blastomere number

## Blastomere size

4 = regular, even division  
 3 = <20% difference (blast diam)  
 2 = 20-50% difference  
 1 = >50% difference  
*after Hardarson et al 2001<sup>33</sup>*

## Fragmentation

4 = <10% frags by volume  
 3 = 10-20%  
 2 = 20-50%  
 1 = >50%  
*after van Royen et al 2003<sup>35</sup>*




## Elective Single Embryo Transfer: Guidelines for Practice British Fertility Society and Association of Clinical Embryologists

RACHEL CUTTING<sup>1</sup>, DAVE MORROLL<sup>2</sup>, STEPHEN A. ROBERTS<sup>3</sup>, SUSAN PICKERING<sup>4</sup>, & ANTHONY RUTHERFORD<sup>2</sup> (ON BEHALF OF THE BFS AND ACE)

<sup>1</sup>Centre for Reproductive Medicine and Fertility, Jessops Wing, Sheffield Teaching Hospitals NHS Foundation Trust, Sheffield, UK, <sup>2</sup>Reproductive Medicine Unit, Leeds Teaching Hospitals NHS Trust, Clarendon Wing, Leeds General Infirmary, Leeds, UK, <sup>3</sup>Health Methodology Research Group, University of Manchester, Manchester, UK, and <sup>4</sup>Edinburgh Fertility & Endocrine Centre, Royal Infirmary of Edinburgh, Edinburgh, UK.

# Blastocyst grading

<b>Expansion Status</b>	<p>1 = Early blastocyst; blastocoel less than half the volume of the embryo, little or no expansion in overall size, zona pellucida (ZP) still thick</p> <p>2 = Blastocyst; blastocoel more than half the volume of the embryo, some expansion in overall size, ZP beginning to thin</p> <p>3 = Full blastocyst; blastocoel completely fills the embryo.</p> <p>4 = Expanded blastocyst; blastocoel volume now larger than that of the early embryo. ZP very thin</p> <p>5= Hatching blastocyst; trophectoderm has started to herniated through the ZP</p> <p>6 = Hatched blastocyst; the blastocyst has evacuated the ZP</p>
<b>ICM grading</b>	<p>A = ICM prominent, easily discernible and consisting of many cells, cells compacted and tightly adhered together</p> <p>B= Cells less compacted so larger in size, cells loosely adhered together, some individual cells may be visible</p> <p>C = Very few cells visible, either compacted or loose, may be difficult to completely distinguish from trophectoderm</p> <p>D = Cells of the ICM appear degenerate or necrotic</p> <p>E = No ICM cells discernible in any focal plane</p>
<b>Trophectoderm</b>	<p>a = Many small identical cells forming a continuous trophectoderm layer</p> <p>b = Fewer, larger cells, may not form a completely continuous layer</p> <p>c= Sparse cells, may be very large, very flat or appear degenerate</p>



# Which cryoprotectant?

- no clear evidence to favour any one system
- optimise system within each laboratory
- anecdotally, DMSO favoured for oocytes?

# How and when to control temperature?



- optimised at RT°C (constant in all seasons)
- permeability of CPAs into cells is temperature-dependent
- rate of ingress of CPA can impact on outcomes
- warm solutions to RT °C at least 30 minutes before use
- switch off heated stages or use ambient plate
- do not leave over light source

# Are media volumes important?

- yes ... and no!
- people use small volumes successfully but...
- control of osmolarity changes due to evaporation
  - dispense only immediately before use
  - be aware of effects of airflow in workstation
- pipetting
- number of oocytes/embryos handled at one time

\*\*\* most new users improved outcomes using larger volumes \*\*\*


# Vitrification Cooling Protocol: ORIGIO






# Vitrification Cooling Protocol: SAGE

Fig. 1



**KEY:**  
 ES = Equilibration Solution  
 VS = Vitrification Solution  
 ➡ = transfer embryo to next drop  
 T = Top of drop  
 C = Center of drop  
 B = Bottom of drop


Fig. 2



**KEY:**  
 ES = Equilibration Solution  
 VS = Vitrification Solution  
 ➡ = transfer embryo to next well



# How long should we equilibrate?



L. Parmegiani - 2011


- stage-dependent
- 5 – 15 minutes
- options
  - observe and wait to see full (90%) re-expansion
  - establish median (fixed) time for your laboratory
- blastocysts – effect of collapse
  - collapsed: use 5 minutes
  - non-collapsed: check for re-expansion of cells NOT blastocoel

# Should we collapse blastocysts?

Reproductive Biology and Endocrinology 2009, 7:99


<http://www.rbej.com/content/7/1/99>

- some clinics doing well without collapsing
- preferred options
  - puncture with ICSI needle
  - laser
  - small tip (micropipetting)
- generally advised
- wait for 50% shrinkage and then move straight to VM



**Figure 1**  
**Different pre-vitification interventions for blastocysts.** A. Assisted hatching: An opening is created in the zona using laser pulse B. Needle blastocoele puncture: A needle is passed through the zona and blastocoele and retracted allowing the blastocoele fluid to freely leak. C. Laser blastocoele puncture: laser pulse creates an opening in the zona and a small defect in the trophectoderm causing the blastocoele to leak. D. Blastocoele aspiration: An injection needle is introduced into the blastocoele and blastocoele volume is sucked out. E. Micropipetting: Passing the blastocysts through a narrow pipette would crack the zona and compress the blastocoele to leak through the cracked zona.


# Should we collapse blastocysts?



video courtesy of  
Lyndon Miles, Cardiff

**CRGW**  
Centre for Reproduction and Gynaecology Wales

# Should we collapse blastocysts?




video courtesy of  
Lyndon Miles, Cardiff

**CRGW**  
Centre for Reproduction and Gynaecology Wales


# Should we collapse blastocysts?

Table 2. Effect of blastocoele aspiration and spontaneous hatching on expanded blastocysts postwarming DNA integrity.

Parameter	Expanded blastocysts			
	Fresh	Vitrified		
	Control (n = 14)	No intervention (n = 22)	Blastocoele aspirated (n = 22)	Spontaneously hatched (n = 12)
Total blastomeres mean $\pm$ SD	48.00 $\pm$ 15.66	46.91 $\pm$ 15.05	39.95 $\pm$ 12.84	62.00 $\pm$ 13.10
Survival	100%	90.9%	100%	100%
DNA integrity index mean $\pm$ SD	95.47 $\pm$ 4.25	77.61 $\pm$ 16.65 ( $P < .001$ ) <sup>a</sup>	90.08 $\pm$ 5.59 ( $P = .003$ ) <sup>a</sup> ( $P = .001$ ) <sup>b</sup>	88.45 $\pm$ 4.35 ( $P = .001$ ) <sup>a</sup> ( $P < .012$ ) <sup>c</sup> ( $P = .35$ ) <sup>d</sup>




# Vitrification step



5 secs



30 secs




60 secs

- do not pre-load pipette with VM/VS
- ensure complete mixing
- care with timing before loading

# Choice of carrier: open or closed?







- loading
  - use a holder
  - minimal volume
  - numbers
  - maximal cooling
  - timing (not too soon so evaporation significant)






# Vitrification







# Storage conditions


- liquid or vapour phase
- care with handling
  - movement to small vessels prior to warming
  - audits
- even brief exposure to temperature above glass transition temperature could be disastrous



# Warming



# Vitrification Warming Protocol: Origio



warming step critical – care with volumes; otherwise use smaller drops ( $\pm$  oil)



# Vitrification Warming Protocol: Sage



Fig. 1







Fig. 2



**KEY:**  
**1M WS** = 1 M Sucrose Warming Solution ➡ = transfer embryo to next drop  
**0.5M WS** = 0.5 M Sucrose Warming Solution **T** = Top of drop  
**MS** = MOPS Solution **B** = Bottom of drop

**KEY:**  
**1M WS** = 1 M Sucrose Warming Solution ➡ = transfer embryo to next well  
**0.5M WS** = 0.5 M Sucrose Warming Solution  
**MS** = MOPS Solution

# Vitrification Warming Protocol



Published in final edited form as:

*Cryobiology*. 2009 August ; 59(1): 75–82. doi:10.1016/j.cryobiol.2009.04.012.

## The Dominance of Warming Rate Over Cooling Rate in the Survival of Mouse Oocytes Subjected to a Vitrification Procedure\*

Shinsuke Seki and Peter Mazur\*

Fundamental and Applied Cryobiology Group, Department of Biochemistry and Cellular and Molecular Biology, The University of Tennessee, Knoxville, Tennessee 37932-2575

# Vitrification Warming Procedure

- prepare in advance
- first step must be at 37°C
  - keep warming medium at 37°C (volume; oil/water in outer well)
  - shortest possible time from liq N2 to WM



- remaining steps at RT°C or 35-37°C depending on kit used

# Comparison of Slow Freezing & Vitrification

**Table 1** Oocyte key performance indicator values.

KPI		Competence		Benchmark
O1	Morphological survival	Freezing	≥50%	75%
		Vitrification	70%	85% (95% for donors <30 years)
O2	Fertilization rate	No more than 10% (absolute; i.e. 10 percentage points) lower than that for the comparable population of fresh oocytes at the centre		
O3	Embryo development rate	Freezing	No more than 10–30% (relative) lower than that for the comparable population of fresh embryos at the centre	The same as for the comparable population of fresh embryos at the centre
		Vitrification	The same as for the comparable population of fresh embryos at the centre	
O4	Implantation rate	No more than 10–30% (relative) lower than that for the comparable population of fresh embryos at the centre		

## ARTICLE

**The Alpha consensus meeting on cryopreservation key performance indicators and benchmarks: proceedings of an expert meeting**

Reproductive BioMedicine Online (2012) **25**, 146–167



# Comparison of Slow Freezing & Vitrification

**Table 3** Embryo key performance indicator values.

KPI		Competence		Benchmark
E1	Morphological survival: fully intact	Freezing	40%	55%
		Vitrification	70%	85%
E2	Morphological survival: $\geq 50\%$ intact	Freezing	60%	85%
		Vitrification	85%	95%
E3	Post-thaw development (including implantation rate) for fully intact embryos	$\leq 10\%$ (relative) lower than that for the comparable population of fresh embryos at the centre		The same as for the comparable population of fresh embryos at the centre

The KPI values are calculated as the proportion of all thawed/warmed embryos.

## ARTICLE

**The Alpha consensus meeting on cryopreservation key performance indicators and benchmarks: proceedings of an expert meeting**

Reproductive BioMedicine Online (2012) **25**, 146–167

# Comparison of Slow Freezing & Vitrification

**Table 4** Blastocyst key performance indicator values.

KPI		Competence		Benchmark
B1	Survival rate	Freezing	70%	85%
		Vitrification	80%	95%
B2	Transfer rate	Freezing	70%	85%
		Vitrification	80%	95%
B3	Implantation rate	≤10% (relative) lower than that for the comparable population of fresh embryos at the centre		The same as for the comparable population of fresh embryos at the centre

The KPI values are calculated as the proportion of all thawed/warmed blastocysts.

[ARTICLE](#)

**The Alpha consensus meeting on cryopreservation key performance indicators and benchmarks: proceedings of an expert meeting**

**Reproductive BioMedicine Online (2012) 25, 146–167**

# How to start vitrification successfully (or optimise your results)

- select only good quality gametes and embryos
- be fastidious with control of timings, temperature, osmolarity and pipetting
- evaluate time for shrinkage and re-expansion in EM (use median time as fixed time)
- load in very small volumes
- take care with handling and storage
- remember that warming is as critical as vitrification steps

